	<b>Title: Website Catalog &amp; specification Document</b>	<b>Document no: NASAN-D&amp;D-D-007</b>
		<b>Rev no: 0</b>
		<b>Date: 22-04-2025</b>
		<b>Details of change: new format</b>

**M-Cardia** is a resting ECG machine. It acquires ECG and prints it on A4 paper. It has e-mail option to send ECG report.

## Technical Specifications

### 1. Performance –

Display : Android mobile\_android 13 and 14

Language : English

Printout : on A4 Paper

Email : Yes

### 2. Power supply

USB powered

### 3. Dimension and Weight


Dimension (in mm) (L x B x H) : 105X65X40 mm

Weight (in Kg) : 150gm

### 4. Environmental

Operating Temp : -5 °C to +50 °C

Operating & Storage : up to 95 % non-condensing

	<b>Title: Website Catalog &amp; specification Document</b>	<b>Document no: NASAN-D&amp;D-D-007</b>
		<b>Rev no: 0</b>
		<b>Date: 22-04-2025</b>
		<b>Details of change: new format</b>

Humidity

## 5. ECG

Frequency response : 0.05 to 100 Hz with notch at 50 Hz

Input impedance : 10 M ohm

Time constant : > 3.2 sec

A/D conversion : 12 bit

Noise : < 20  $\mu\text{V}_{\text{p-p}}$

Leakage current : < 10  $\mu\text{A}$

CMRR : > 120 dB

Patient isolation : Optical isolation.

Lead selection : standard 12 leads

Gain : 5, 10 & 20 mm/mV

Sweep speed : 25.0, 50.0 mm/S.

## Detailed Feature

### ➤ Main screen display

- **Patient** –To go to patient search page
  - Search old patients
  - Register new patient
  - Edit medical history
  - Edit patient detail
  - Link emergency
  - Conduct ECG test

	<b>Title: Website Catalog &amp; specification Document</b>	<b>Document no: NASAN-D&amp;D-D-007</b>
		<b>Rev no: 0</b>
		<b>Date: 22-04-2025</b>
		<b>Details of change: new format</b>

- View log of patient
- **Review tests** – To view list of all conducted tests.
  - view report or check the status of the patient with date and time of the test
  - Facility to Print A4
  - Facility to send email
  - Preview ECG
- **Emergency** –To take emergency ECG
- **Configuration**
- **Details of Configuration**
  - Screen configuration (Default screen)
    - ECG screen
    - Median screen
  - Report Configuration
    - Resting ECG with Measurement
    - Resting ECG with Measurement and Interpretation
  - Short details check box
  - Grid/Graph/No grid option
  - Color, Borders check box
  - Long lead selection from dropdown for A4 report
  - Select emergencies number from dropdown
  - Select sleep time from dropdown
  - Facility to enter hospital name
  - Facility to enter send to mail ID
  - Facility to enter doctor name
  - Cancel & OK button
- **Acquisition display (ECG screen)–From emergency & new ECG button**
  - Simultaneous acquisition of 12 Leads of ECG @ 250 samples per second per channel with ADC resolution of 12-bit.
  - Frequency Response – DSP filter of 0.05 to 150Hz with notch at 50Hz
  - Display of 4 seconds of ECG traces of 12 leads in 6 X 2 with standardization pulse
  - Acquisition screen button
    - Exit button
    - Gain button
    - speed button
    - QRS lead button
    - Save button
  - **Gain button-** allows you to choose desired recording sensitivity (5, 10, 20mm/mv).
    - the choice of changing the gain of the Electrocardiogram from acquisition

	<b>Title: Website Catalog &amp; specification Document</b>	<b>Document no: NASAN-D&amp;D-D-007</b>
		<b>Rev no: 0</b>
		<b>Date: 22-04-2025</b>
		<b>Details of change: new format</b>

- **Speed button-** asaan allows you to choose desired recording sensitivity (12.5, 25, 50mm/s).
  - allows the choice of changing the sweep speed of the Electrocardiogram from acquisition
- **QRS lead button-**
  - To change QRS lead and change lead displayed in magenta color
- **Save button-**
  - This button will get enable after 10 sec.
- **Exit button-**
  - Exit from the acquisition screen
- **Median Screen display-**
  - **Zoom selected median**
    1. Facility to change P+ start, P+ end, P – end, Q start, R Start, R end, S end, r' end, s' end, T Start, T End, U start, U end points.
  - **Displays following Measurement**
    1. Amplitude (in V)
    2. Duration of P,Q,R ,S,T in ms
    3. Interval of PQ,PR,QRS,QTC in ms
    4. Axis of P,QRS and T measurement in degree
    5. Heart rate
  - **Displays following ECG interpretation**
    1. Rhythm statements
      - a. Sinus Rhythm / Junctional or NonSinus Rhythm suspected / Sinus Bradycardia / Sinus Tachycardia
    2. PR interval statements
      - b. PR is normal / PR is long / PR is short
    3. QRS interval statements
      - c. Normal QRS Width / Wide QRS
    4. QT interval statements
      - d. Normal QT interval / QT Interval is prolonged / QT Interval - Little clinical significance
    5. QRS Axis statements
      - e. QRS Axis is normal / Right axis deviation / Left axis deviation / QRS Axis is indeterminate
    6. ST abnormalities statements
      - f. Normal ST level / ST Elevation / ST Depression
    7. T wave inversion
      - g. Lead names of leads showing T wave inversion
  - **Edit Measurement & Interpretation facility**
  - **Facility to Enter comment**

\*Due to our continuous product development program, the specification may change without prior notice.