


Pages 1 / 1	Specifications Product: AMRUT	
-----------------------	--	---

1. General

Dimensions : 7.5cm x 20cm x 24cm

Weight : Approx.2.2 Kgs.

Controls : On/Off button ,Red shock button

1. Battery

Type :15V/4.2 AH,lithium manganese dioxide, disposable long life primary battery

Capacity : minimum 200 shocks of 200 Joule or 4 hours

Standby life : Four years typical when battery is installed by the Install by date

2. Environmental

Operating Temp : -5°C to +50 °C

Humidity : Up to 95 % non-condensing

3. AED

Wave form : Biphasic truncated exponential

AED Shock Series : Energy level 200J

Charge Command : Automatic after identifying shockable rhythm

Shock Command : Front panel button, shock

Defibrillation Electrode : Self adhesive, adult pads

Voice and Text prompts : Extensive text and audible messages to guide user

Shock to Shock Cycle time : Typically less than 20 seconds between shocks in a series.

* Due to our continuous product development program, the specification may change without prior notice.