



AMRUT First Indian AED

Key Features

- ✓ Detailed voice prompts and visual icons for guidance to user
- ✓ Battery capacity : 200 shocks or 4 hours
- ✓ 20 sec delay time between two shocks
- ✓ Evaluates patient ECG to determine if rhythm is shockable

For Technical Specification refer :

 <http://www.nasanmedical.com/pdf/specs/Amrut-S.pdf>

NASAN MEDICAL ELECTRONICS PVT.LTD

Plot No. 147, "Rangali" Society, Opp. Karve Statue, Kothrud, Pune - 411 038.

Phone : +91 82370 32389 E-mail : sales@nasanmedical.com Website : www.nasanmedical.com