

Pages 1 / 1	Specifications Product: Sanjeevani 1001	
-----------------------	--	---

1. General

Display : 7" colour TFT
Display area : 150 mm x 90 mm
Dimensions : H240 x W 230 x D200 mm (Max) without Paddle
Weight : Approx 5.0Kgs

2. Power Supply

AC mains through adaptor
Input : 150 - 270 VAC, 50 Hz, 1 Ø
Output : 17 @ 4 A.
Battery : 11.1 V 4Ah Lithium-ion .Rechargeable
Back-up time : 3 hours min. on full charge,
70-100 shock of 150 Joules deliveries (Full charge).
Power consumption : 25 Watts (max.) in ECG mode

3. Environmental

Operating Temp. : -5°C to +50 °C
Humidity : up to 95 % non-condensing

4. Defib

Operating Mode : Manual
Waveform : Biphasic
Energy : 2 to 200 joules in 12 steps selection
Synchronous Cardio Version : Energy delivery begins within 60 ms of the QRS Peak
Charging Time : Less than 8 seconds
Paddle Assembly : External paddle assembly with coil retractable high voltage cable. Adults & Pediatric paddle integrated to same handle
Check Mode : Manual test mode
Operation : By front panel keyboard and through Encoder, Unit ON-OFF by front panel push button

* Due to our continuous product improvement programme, features can be enhanced.