

# Biphasic Defibrillator



Optional  
AED  
Pacer  
SpO2  
NIBP  
EtCO2  
Internal Paddle

\*Also available in 360 Joules

## SANJEEVANI 1006

### Key Features

- ✓ 7" Colour Display with built in printer
- ✓ Charging time 8 sec on mains
- ✓ Offers variable energy level selection (2,3,5,10,20,30,50,70,100,150,200 Joules)
- ✓ Battery capacity 70-100 shocks on single charge
- ✓ Adult & pediatric paddles support
- ✓ Compact & Light weight
- ✓ Print out on 2" thermal paper
- ✓ Voice prompt facility


#### Premium Defibrillator Model

- Charging time less than 4 second
- SD card storage facility
- 250 shock events with 10 second of ECG data
- Contact Indicator

For Technical Specification refer :

[http://](http://www.nasanmedical.com/downloads/specifications/Specification-Sanjeevani-1006.pdf)

[www.nasanmedical.com/downloads/specifications/Specification-Sanjeevani-1006.pdf](http://www.nasanmedical.com/downloads/specifications/Specification-Sanjeevani-1006.pdf)



\*Automatic Print out generated after every shock delivery

**NASAN MEDICAL ELECTRONICS PVT. LTD.**

C wing, Shubham Heights Off Mumbai - Bangalore Highway Service road, Warje Pune - 411058

Phone: +91-8237032389

Email: [sales@nasanmedical.com](mailto:sales@nasanmedical.com)

[www.nasanmedical.com](http://www.nasanmedical.com)