



- **OPTIONAL**
- AED
- Pacer
- SpO2
- NIBP
- EtCO2
- Internal Paddle

*Also available in 360 Joules

Sanjeevani 1006


Key Features

- ✓ 7" Colour Display with built in printer
- ✓ Charging time 6 sec on mains
- ✓ Offers variable energy level selection (2,3,5,10,20,30,50,70,100,150,200 Joules)
- ✓ Battery capacity 70-100 shocks on single charge
- ✓ Adult & pediatric paddles support
- ✓ Compact & Light weight
- ✓ Print out on 2" thermal paper
- ✓ Voice prompt facility

Technical Specifications

Premium Defibrillator Model

- Charging time less than 4 second
- SD card storage facility
- 250 shock events with 10 second of ECG data
- Contact Indicator



*Automatic Print out generated after every shock delivery



[Click Here For Technical Specification](#)



*Due to our continuous product development program, the specifications may change without prior notice.

NASAN MEDICAL ELECTRONICS PVT. LTD.

C Wing - Shubham Heights, Off Mumbai - Bangalore Highway Service Road, Warje, Pune - 411058

Phone : +91 8237032389 | E-mail : sales@nasanmedical.com

www.nasanmedical.com